



The Czech Academy  
of Sciences

# Scanning Enhancement of STM-tungsten Probes by Applying Colloidal Graphite Coatings

Mohammad M. Allaham, M.Sc.


Central European Institute of Technology, Brno University of Technology, Purkyňova 123, 612 00 Brno, Czech Republic.

Group of e-beam lithography, Dept. of Electron and plasma technologies, Institute of Scientific Instruments of CAS, Královopolská 147, 612 64 Brno, Czech Republic.

# Introduction

Scanning tunneling microscopy (STM) is an atomic level high-resolution surface analytical technique.

STM is a vacuum tunneling technique that uses a good conductive metallic sharp tip (tungsten, platinum/iridium, iridium, ...) to scan the topography of a conductive sample surface.




STM of Graphene

Tungsten nano cold field emission tips is one of the most used probes in STM.


However, tungsten STM-probes have some limitations:

- 1- Requires ultra-high vacuum conditions.
- 2- Requires in-situ cleaning before a scanning process.
- 3- Sensitive against adsorbates and is quickly oxidized.




Tungsten cold field emission STM probe.

# Introduction




First approximation for the tunneling Schottky-Nordheim potential energy barrier between the sample and probe

# STM from uncoated W cathodes in SEM chamber



Result of the STM scan attempt using a tungsten cold field emission tip ( $P = 10^{-2}$  Pa).

# Limited surface conductivity of W



XPS spectra of cleaned tungsten after different Ar500+ ion etching steps

# Limited surface conductivity of W


—Pre-etching —5 keV/2 min —10 keV/2 min —20 keV/2 min

— Pre etching  $\phi_1 = 3.65 \text{ eV}$

— 5 keV/2 min  $\phi_2 = 3.78 \text{ eV}$

— 10 keV/2 min  $\phi_3 = 4.03 \text{ eV}$


— 20 keV/2 min  $\phi_4 = 4.73 \text{ eV}$



UPS spectra of cleaned tungsten after different Ar500+ ion etching steps




# Why colloidal graphite?



SEM of colloidal graphite surface deposited on glass


# Why colloidal graphite?



Raman spectroscopy for a graphite sample surface compared to amorphous polymer graphite




# Why colloidal graphite?




Raman spectroscopy peaks analysis for polymer graphite (a and b) and for colloidal graphite (c and d)

# Why colloidal graphite?








Reflected electron energy loss for colloidal and polymer graphite along with HOPG

# Why colloidal graphite?



UPS spectrum of colloidal graphite

# Samples prepared





# Sample prepared





Fig.10. Tungsten-graphite cathode



SEM calibration sample





Fig.9. Glass-graphite cathode




# Setup – FEM




Zero series resistor field emission microscope

# Setup – STM




[www.nenovision.com](http://www.nenovision.com)




Nenovision litescope AFM in SEM

# Field emission characteristics of glass-graphite cathodes



The field emission characteristics of the glass-graphite microtips.

# Field emission characteristics of glass-graphite cathodes




The field emission current-time characteristics of the glass-graphite microtips at pressure of 1 Pa.




# STM with a W-Graphite cathodes

(c) STM-top view




(d) STM-bottom view



STM measurement with a W-Graphite cathode ( $P = 10^{-2}$  Pa)




# STM with a W-Graphite cathodes




# Summary

(a)  
No applied electrostatic field



(b) With applied electrostatic field



## Summary

1. Tungsten cold field emission STM probes have limited performance outside UHV environments.
2. XPS and UPS showed a critical requirement for in-situ cleaning procedure in FEM and STM instruments.
3. Colloidal graphite showed good performance as coating material for cold field emission applications.
4. W-Graphite showed good cold field emission performance.
5. W-Graphite cold field emission STM probes showed good STM performance inside SEM chambers for topography studies.
6. Graphenic structure of graphene related materials is believed to be the reason behind the linear behavior in the corresponding current-voltage characteristics.

# Thank you for your attention.

## Questions?

**Mohammad M. Allaham, M.Sc.**

Ph.D. candidate  
Electron and Plasma  
Technologies



**ISI**  
**CAS**

Institute of Scientific  
Instruments  
The Czech Academy  
of Sciences

M: +420 735 039 950  
E: [allaham@isibrno.cz](mailto:allaham@isibrno.cz)  
[www.isibrno.cz](http://www.isibrno.cz)

Královopolská 147  
612 64 Brno  
Czech Republic